

Polierer Polisher

Silikon silicone












Silikonpolierer sind ideal zum Abtragen, Glätten und Polieren von Acryl und Prothesenkunststoffen. Ebenfalls geeignet sind sie für die schonende Nagelpflege.

In der Bellflex-Ausführung eignen sie sich für einen kühlen Schliff und eine effektive Politur von Silber, Platin und Titan sowie für dentale Korrekturen ohne Glanzbrand. In der Belldia-Ausführung sind sie ideale Diamantpolierer für Hochglanz auf Keramik und Metall.

Silicone polisher are ideal for removing, smoothing and polishing acrylic and overdenture plastics. Likewise suitable for gentle nail care.




















In the Bellflex version suitable for a cool finish and effective polish of silver, platinum and titanium as well as for dental corrections without glaze firing.

In the Belldia variant they are ideal diamond polishers for providing a high gloss on ceramics and metal.

Art.-Nr. unmontiert art.-no. unmounted	Art.-Nr. montiert art.-no. mounted	Härte grade	Schaft shank Ø mm	 mm	 mm	Mandrell für unmontiert auf Anfrage Mandrel for unmounted on enquiry	 10	 100
Silikonpolierer für Gold- und Silberlegierung silicone polisher for gold and silver alloy						8.000/min		
SP193-01		grob coarse		22	6			
SP193-02		grob coarse		22	3			
SP193-03		grob coarse		17	3			
SP193-04		grob coarse		22	4			
SP193-05		grob coarse		18	3,5			
SP193-06		grob coarse		6	24			
SP193-07		grob coarse		7	20			
	SP193-08	grob coarse	2,34	14,5	2,5			
	SP193-09	grob coarse	2,34	14,5	5			
	SP193-10	grob coarse	2,34	14,5	2,5			

Polierer Polisher










Silikon silicone






Art.-Nr. montiert art.-no. mounted	Härte grade	Schaft shank Ø mm	 mm	 mm	
Silikonpolierer für Prothesenkunststoffe silicone polisher for prosthesis plastics					 10 100 VE Stück/PU pieces
					 6.000/min
S195-1G	grob coarse	2,34	12	22	
S195-2G	grob coarse	2,34	10	20	
S195-3G	grob coarse	2,34	5,6	16	
S195-4G	grob coarse	2,34	10	24	
S195-5G	grob coarse	2,34	15	17	
S195-1M	mittel medium	2,34	12	22	
S195-2M	mittel medium	2,34	10	20	
S195-3M	mittel medium	2,34	5,6	16	
S195-4M	mittel medium	2,34	10	24	
S195-5M	mittel medium	2,34	15	17	
S195-1F	fein fine	2,34	12	22	
S195-2F	fein fine	2,34	10	20	
S195-3F	fein fine	2,34	5,6	16	
S195-4F	fein fine	2,34	10	24	
S195-5F	fein fine	2,34	15	17	

Polierer Polisher

Bellflex / Belldia-Diamant

Bellflex / Belldia diamond

Art.-Nr. unmontiert art.-no. unmounted	Härte grade	Farbe color	 mm	 mm	Mandrell für unmontiert auf Anfrage Mandrel for unmounted on enquiry	 10  100 VE Stück/PU pieces
Bellflex_effektive Politur auf Titan, Platin und andere Metalle effective polish on titan, platin and other metal  10.000/min						
BF209-G	grob coarse		22	3		
BF209-M	mittel medium		22	3		
BF209-F	fein fine		22	3		

Art.-Nr. unmontiert art.-no. unmounted	Art.-Nr. montiert art.-no. mounted	Härte grade	Schaft shank Ø mm	 mm	 mm	
Belldia_Hochglanz auf Keramik (auch für Metalle geeignet) for high gloss on ceramics, also on metal  10.000/min						
BD209-1F		fein fine		17	2,5	
	BD209-2F	fein fine	2,34	11	2	
	BD209-3F	fein fine	2,34	4	13	
BD209-1M		mittel medium		17	2,5	
	BD209-2M	mittel medium	2,34	11	2	
	BD209-3M	mittel medium	2,34	4	13	